

# SKIN COOL

Hot & Cold & EMS rejuvenation instrument




Thank you for buying this product  
Please read the instruction carefully before using.  
Thank you for your cooperation!

# Contents

Principle-----	1
Host diagram -----	2
instructions of a course of treatment-----	2
Function parsing -----	4
Treatment effect-----	6
Operation diagram of the face -----	7
Notice-----	7
Service and repair-----	9

SKIN COOL has three features. Mainly used as "DNA cells injection course of" Korean medical beauty. Now in more than 700 thousand have been used SKIN COOL instrument, can use the COOL operation, at the same time inject the essence, make the skin cells to absorb the products up to 80%, the penetration ratio is more than 20 times, SKIN COOL instrument has 10 effects:

1. Promote lymph circulation.
2. Minus the expression lines, light bitter lines.
3. Remove the dark yellow skin and edema.
4. Minus double chin.
5. Effectively eliminate dark circles and eye pouch.
6. Shrink pores, promote blood circulation.
7. Anti wrinkle, enhance eyes round.
8. Activate collagen layer to promote skin cell regeneration.
9. Enhance more skin elastic, more compact.
10. Enhance facial contours, promote fat dissolution.



## SKIN COOL indicator diagram

### EMS + COOL & HOT course:

HOT+EMS: HOT import + EMS can receive more than **Course operation**

### tips:

1.40000 Dalton solution, through the temporary microporous membrane directly absorbed by the skin cells.

2.COOL+EMS:Use of - 10 ° freezing function, make the skin cells to absorb and lock and after shut down.

Note: it's best to do HOT thermal + EMS electroporation, redo frozen injection + electroporation, so it can make moisture into the skin nutrition is well locked and contractive pore effect.

By injecting -10 degrees of freezing, the principle and essence to be worn out pulsation frequency cells into cells, promote products increased 80% utility, to enhance the repair and improve skin problems.

## **HOT injection technique to cell** —using serum or essence injection

Cell injection efficacy, reduce the fine lines, promote cell regeneration, increase skin absorption capacity;

Penetration of nutrition essence, stimulate cell vitality, compact skin relaxation;

Nourish cells, increase blood circulation, stimulate metabolism.

## **COOL and extremely slow injection technology**

Anti allergy, relieve, close pores, restore skin, sedative effect,

Strengthen skin elasticity, improve skin regeneration capacity.

Increase the absorption capacity of the skin, keep the skin moist and soft.

## **Shrink pores and bright nose method**

Collagen hyperplasia, increase elasticity. Coarse pores, such as the nose next to the face position, has obvious effect in the course of contraction. The probe close to the nose, forced to massage as a circle to the nose, the action is repeated 3 times.

## **Invincible V face method**

Tighten facial muscles, remove edema, compact skin, shape V face, while accelerating the blood circulation, to remove yellow skin, skin whiten and turn red, make pore more compact.

The probe is close to the chin, outward circle to the ear beads along beside the ear pull back to the ear, the ear down to the clavicle, and then in the opposite direction, the clavicle light force pulled up into the jawbone, then pulled forward back to the chin position, the action is repeated 3 times.

2.The probe is close to the mouth, outside a circle to the ear, and then pull back along the ear to ear, from the ear to the clavicle to drop down, and then in the opposite direction, the clavicle light force pulled up into the jawbone, then pulled back to the mouth forward position, the action is repeated 3 times.

3.The probe close to the zygomatic bone (near the nose), outward circle to the ear, and then pull back along the ear to ear, from the ear to the clavicle to drop down, and then in the opposite direction, the clavicle light force pulled up into

the jawbone, then pulled forward back to view the bone, the action is repeated 3 times.

### **Eye rims & Eye pouch**

Make collagen regeneration, remove edema, skin, remove eye pouch, eye wrinkles disappear naturally. At the same time stimulate the blood circulation, reduce dark circles, so that the eyes become beautiful and moving. The probe close to the pouch position, from the eyes of about 1CM, pull to the temple, the action is repeated 3 times.

Note: pouch part without massaging looped, simply pull by the canthus keep along the temple up eye bags and keep 1 cm distance with eyes, avoid to cause not unobstructed.

### **Wrinkles removal and rejuvenation**

No matter how good your maintenance, a deep hard lines, can make you decades older, using SKIN COOL instrument help regenerate muscle layer of collagen, increase elasticity, decrease hard wrinkles, stay young forever.

Methods:

1. Place the probe close to chin, in order to draw the C word method, gently on to philtrum position, place on the philtrum position, for the C word method again, light force landed on the chin, the whole movement repeat 3 times.
2. Place the probe wrinkles location close to the mouth, gently outward looped to the top of the wrinkles beside (nose), light force down again back to the corners of the mouth wrinkle position, the whole movement repeat 3 times.

### **Neck lines double chin**

To eliminate double chin, chin contour clear, how the recent low angle shots are not afraid, and the neck skin becomes more compact, dodge neck lines, thereby blocking your young secret.

Methods: place the probe close to double chin, back behind the ears, and from behind the ears down to the clavicle artery location, and then in the opposite

direction, by the collar bone light force to pull up to the jawbone, then pull back to double chin forward position, the whole movement repeat 3 times.

### **Lymphatic facial surgery**

1.SKIN COOL instrument can dredge lymphoid, enhance the capacity of detoxification and remove edema, to make the future beauty course play the strongest effect.

Methods: place the probe to the clavicle artery location to outgoing circle to behind the ears, and favorable along the aorta down to the starting point (collar bone nearly artery location), the whole movement repeat 3 times.

### **Wrinkles removal and rejuvenation**

Before the massage, must evenly apply to the face and neck with the essence.

Method of use: you can use to go once every other day, 10-15 minutes every time. Suggested two or three times a week, the skin will have a significant effect.

#### **Operation instruction:**

1. Connect the host with the connection lines one by one according to the diagram.
2. turn on the power switch.
3. press start switch.
4. Selection of extreme temperature HOT model, the temperature value is selected as 0-40°C, painted cream or liquid on the face before operation, in accordance with the personal benefit from setting the appropriate temperature.
5. Select the RF EMS output intensity: pressing (press open, then press for recovery) - RF rotary button adjustment strength; strength should be selected according to the individual benefit of determined, for the first time use, suggest moderate or weak intensity, then gradually increase the intensity.
6. After the completion of the HOT operation mode, select the mode of operation of COOL, -10°C to 5°C, the operation method is the same with the HOT mode.
7. After the operation is finished, please close the operation interface and cut off the power supply.

**Product name:** SKIN COOL

**Rated Voltage:** 110V-220V common use

**Out frequency:** 40.68HZ-----50HZ

**Power:** 40W

**Tempature:** Cold -10°, Hot 40°

**Size:** 220 (W)\*220(D)\*80(H)MM

**Weight:** 1.5KG

## Service and repair

1. This machine from the date of purchase (with a note) host a year warranty. Quality problems within one year have free repair, such as the problems after one year will receive maintenance costs. Other has been guarantee deadline, can be responsible for maintenance, but charge spare parts and consumables cost fee. If need supplement or purchase, please contact our company.

2. Causes the following use of personal failure, does not provide free warranty service.

- 1) The unauthorized entry, modification caused the failure of the product;
- 2) Use accidentally blows, falls caused by the fault;
- 3) Lack of reasonable maintenance caused the failure;
- 4) Did not at the correct guidance of instructions caused by the fault.

Maintenance Report		
Product Name : _____ Date of Purchase : _____		
Sub-office: _____ Add: _____ Tel: _____		
Customer Name: _____ Add: _____ Tel: _____		
Maintenance Date	Trouble Reason	Trouble Remark