

RED LIGHT THERAPY PAD

USER MANUAL

Product Specification		Picture
Product Size	Boot length about 1100MM	
Product Power	120W	
LED QTY	528pcs	
LED Type	5050 SMD	
Spetrums	660nm:850nm =1:2	
Irradiation area	0inch>45mW/cm ²	
Adaptor	12V 10A	
Timing	10~90min,9 levels of timing	
Main Materials	Diving material	
Lifespan	50000 Hours	

Introduction:

The red light therapy mat is based on medical red light therapy 660nm and 850nm red infrared light to activate body cells to promote cell metabolism, 660nm facilitate cells collagen growth for beauty care. 850nm infrared light deeply penetrate into the subcutaneous tissue to promote blood circulation, and regulate the nervous system of various parts for healthy care.

Functions:

Red 660nm combined with infrared 850nm promote blood circulation, increase skin collagen, accelerate metabolism, pain relieve, and promote wound healing.

Sports and fitness:

Heal sports injuries and improve cell viability.

Using tips:

1. Inset the DC port into the product connector;
2. Plug in power supply;
3. Power supply: Adapter

Controller Instructions

SMART MULTI-FUNCTION CONTROLLER

Auto-off & Multiple Mode Options

Turn On/Off



P1



P2



P3



P4



P5

Brightness Setting

P1/P2/P3/P4/P5

Pulsed Light

10Hz/40Hz



First Mode

10Hz Pulse Mode



Second Mode


40Hz Pulse Mode



Timing Function

10/20/30/40/50/60/70/80/90Mins


Controller function :


1. Button  +/- adjust the phototherapy with hot compress function:


Red 660nm and near infrared 850nm all brightness.

(Near-infrared light is an invisible wavelength)

2. Button  ON/OFF;

3.Button  control the pulse light 10Hz and 40HZ pulse light adjustable;

4. Button  +/- Energy Increase button: P1/P2/P3/P4/P5;

5. Button  +/- adjust the time: 10/20/30/40/50/60/70/80/90 Minutes.

Notice:

1. If the body feels uncomfortable with the red light therapy, please adjust the low intensity or stop using it immediately;
2. If the skin appears abnormal, please stop using it;
3. Please avoid eyes looking at the red light directly;
4. To ensure consistent optimal performance, always use genuine parts supplied by our company. The use of non-genuine parts may result in product damage, reduced performance, safety hazards, and void the warranty.
5. It is forbidden to soak or wash this product, please turn off the power in time when worked;
6. Do not disassemble or repair without authorization, if there is any quality problem, please contact the seller in time;
7. Pregnant women are forbidden to use;
8. Patients with heart disease are prohibited to use .

Thanks

Disclaimer

1.This product is not intended to diagnose , cure or prevent disease . The company makes no claims, representations or warranties regarding the ability of this product to cure any physical, skin , or mental conditions .A qualified health professional should be consulted with regard to anycondition requiring medical attention.

2.The company makes no legal liability for claims made due to safety issues caused by improper use.

Warranty Card

Product Name			
Product Mode No.			
Date of Production			
Warranty Content			
User Name			
Contact Phone Number			
Mail Address			
Dealer Autograph		User Signature	