

颜色红

C30 M100 Y100 K0

藕荷色

C11 M31 Y33 K0

香槟金

C24 M30 Y50 K0

灰色

CO MO YO K60

材质：128克双铜

说明书为折页

单页尺寸：80\*140mm

双面印刷


7 Colors Photon Skin-Rejuvenation Facial Mask



Read all the instruction before use

Product Part Figure

Note: The infrared light source is not visible light, 30% of the light spots of the inverse age phototherapy beauty instrument are infrared light, Therefore, no visible light can be observed in some light spots after power-on. This is normal. These infrared lights To determine whether the electricity is working properly, heat can be performed by infrared detector or laser thermometer.



1. LED Therapy Mask
2. User Manual
3. USB cable
4. Power Controller

parameter

Product name: 7 colors Photon Skin-Rejuvenation Facial Mask

Model: MR254D

Working Current: 3A

Battery Capacity: 4000mAh

Size: 102.4 \* 28 \* 25.5mm

Charging Time: About 8H

Use Time:120min

Voltage: 3.7V

Power: 12W

Material: PP

Weight: 335g

Color: Diamond White

Product Instructions

When the controller is connected to the USB port state of the mask light:

1. Directly press the switch button, the battery indicator will light up and go out after 2 seconds;
2. Long press the switch button for about 1.5 seconds, the battery indicator light will light up, and the mask light will light up.

Charging or use process

1. When charging, the 4 battery indicator lights flash alternately from left to right; When the battery is fully charged, the 4 battery indicator lights will light up for 1 minute at the same time, and then go out; When the battery is fully charged, unplug the charging head and plug it in again. The battery indicator lights flash from left to right for a few minutes in turn, then the 4 battery indicator lights will light up and go out after 1 minute.

2. When the controller is connected to the mask in use, the 4 battery indicator lights from the right to the left side is off, it means the battery is decreasing; When the battery is low, the leftmost grid light flashes red.

Remarks: When the battery is low, the first light on the left of the battery indicator flashes. This process The mask light is still on at normal brightness, and after 10 minutes, the battery indicator light is off, and the mask light It also goes out. When not charging, press the switch, the first grid light on the left of the battery indicator light flashes It blinks for 5 seconds and then goes out, and the mask light will not light up again during this process. The charging current should be upper than 800mah

Fast mode

Press the button twice quickly to switch to fast mode. In this process, the battery indicator light flashes 3 times from left to right, then the battery indicator is displayed normally, and the mask light turns on in turn

Eight modes:

the light turns on for the first time (blue light)/press again (cyan blue) /press again (white)/press again (green light)/press again (yellow light) /press again (red light) )/Press (purple light) circulate this process have 20 minutes (from the first boot to the end, regardless of whether you press any mode, the time will decrease (for example, when the blue light works for 3 minutes when the blue light is turned on, it will switch to cyan light). The remaining time is only 17 minutes, which is similar to this decreasing) change, and any mode will not re-time) (In the light therapy state, press and hold the power button for about 1.5 seconds, the mask light will go out, and the battery indicator light will go out after 2 seconds ( The controller has no memory function).



Recommended sequence of skincare procedures



1. Please wash your face gently with cleansing products and wipe dry. (Figure 1)
2. Connect the external wiring USB of the inverse age phototherapy beauty instrument to the controller interface, keep the instrument and control Controller connection. (Figure 2)
3. Wear the equipment, long press to turn on the instrument, adjust different gears according to skin needs, enjoy 20 Minutes photon physiotherapy. During the period, you can close your eyes to relax and listen to the music, and it will not affect if you reading and watching movies at the same time. (Figure 3 )
4. After 20 minutes, the power will turn off automatically, remove the instrument, pull out the controller, and close the instrument. (Figure 4 )
5. Continue the daily skin care process, for example: use a regular skin care mask-wash off the mask-use note. (Figure 5 )

Tips

Please stick to it once a night and use it continuously for 7-14 days, you can see obvious results. 28 days is a treatment, three treatment can completely improve the skin base. This product has been clinically verified for long-term use and will prevent skin problems, prevent acne and wrinkle.

Principle and clinical confirmation of phototherapy

- Blue light: Calm the skin, Balanced secretion control grease.
- Cyan light: Metabolism boost energy.
- White light: Decompose stains improve fine lines speed up metabolism.
- Green light: soothes skin, reduces oil secretion and balances skin water oil ratio.
- Yellow light: Decom pose pigment lymphatic drainage improve roughness.
- Red light: Improve skin deep repair skin rejuvenation and wrinkle removal.
- Purple light: Treat acne repair acne marks.

Acne removal example: Figure 6 is before treatment, and Figure 7 is the happening after a course of treatment (28 days) using red and blue infrared mixing the result after the phototherapy treatment, Figure 8 is a photo afer two blue light treatments. Visible blue light high therapeutic effect of physiotherapy on acne.




Anti-aging example: Figure 9 is befoere treatment, Figure 10 is a continuous course of treatment (28 days) daliy use the effect after red light and infrared physiotherapy. It can be seen that the fine lines of the eyes are greatly reduced, the skin is firm, and the skin lower collagen surplus. The combination of red light and infrared can effectively resist skin aging.




Acne effect diagram

Figure 11 is the state before treatment, and Figure 12 is the effect after a blue light physiotherapy cour Figure 13 is the effect after two blue light treatments. It can be seen that blue light has a high curative effect on acne.



Anti-sun old effect picture

Figure 14 is the state before treatment, Figure 15 is the effect of continuous use of red light infrared physiotherapy, Figure 16 is the effect after using red light infrared physiotherapy for 30 consecutive days. Visible fine lines Greatly reduces the firmness of the skin and fills the collagen under the skin. The combination of red light and infrared can effectively resist skin aging.



In February 2004 the famous blu-ray surgery magazine published the blue light for acne bed treatment results. Clinical trials have proved that more than 85% of patients after eight blue light treatments, Facial acne has improved significantly. Even after 3 months after stopping blue light therapy, more than 80% of patients will not relapse. Blu-ray physiotherapy shows unprecedented treatment of acne skin the recovery rate.

In 2005 a paper published in the Journal of laser physiotherapy and beauty pointed out that the combination of red light and infrared light can fundamentally improve facial wrinkles and reduce stains within 12 weeks. Effective clinical treatment up to 81%.

7 Colors Photon Skin-Rejuvenation Facial Mask

