PIAOSY

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PICOSECOND PEN

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Please read the instructions carefully before use and keep the instructions for future reference

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Press the power button for three seconds, the machine starts, the frequency shows 0, and the energy shows 0.

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Press the energy adjustment button and the frequency adjustment button to the appropriate gear.

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Contact the picosecond pen vertically with the skin surface, and align the launching hole with the position to be operated.

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When you need to stop the operation, just release the transmit button.

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Protective glasses must be worn during operation with blue laser.

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ENERGY AND FREQUENCY

There are nine kinds of frequencies and energy in the pen. When the frequency is 1, the number of times of emission per unit time is the least, when the frequency is 9, the number of times of emission per unit time is the most, when the energy is 1, the energy of emission per unit time is the least, and when the energy is 9, the energy of emission per unit time is the largest.

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Avoid consuming ginger which can cause wrinkled scar that is hard to be treated.

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Attention shall be paid to cleaning to avoid exposure and water.

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The instrument cannot be used during the charging process. In case of abnormal heating, the power shall be cut off immediately

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The instrument shall be kept away from the place where children are not easy to touch, and children are not allowed to play.

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It is strictly forbidden to aim the instrument at the eyes, or observe the eyes in case of failure.

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DEAR USER:

We sincerely thank you for purchasing and using our products. We solemnly promise a one-year warranty for To ensure that you can fully share the perfect after-sales service provided by our company, please read this warranty card in detail, Fill in the following information together with the seller and keep it properly.

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Instructions for use of quality assurance card

1. Please fill in the product warranty card accurately;
2. This warranty card is saved by the user;
3. If you need warranty, please show the product warranty card and purchase receipt;
4. This warranty card must be signed or sealed by the seller before it is valid;
5. If any quality problem is found within three days after purchasing the product, ask the dealer for free replacement;

The main engine can be repaired free of charge within one year; the manufacturer is responsible for the repair free of labor cost within ten years (only the original is recovered Material cost).

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