



Switch with display screen



-  **On/Off**
-  **Three modes:**  
660+850 / 660 / 850
-  **Dimming key:**  
5 gear dimming
-  **Timing key:**  
any adjustment from 5 to 30 minutes

## FLEXIBLE & ADJUSTABLE



Portable and Foldable

360°



▶ Rotatable Gooseneck



▶ Adjustable Stand



▶ Stable Base




Height Range From 35" to 67"



When the tops of the stand's three legs are 25.6 inches apart, the stand will be more stable and effectively prevent tipping.

## Verbesserter Drei-in-Eins-Chip


180PCS 850nm Nah-Infrarot-Licht  
360PCS 660nm Rotlicht



## Wie die Rotlichttherapie Funktioniert

Durch die besondere Gestaltung des optischen Systems und High-Power-LED-Lichtquelle, die Strahlungsleistung eines bestimmten Spektrums und Energiedichte der schmalen Frequenz Niedrigenergie-Rotlicht, den menschlichen Körper zu produzieren photochemische Effekte, anstatt thermische Effekte.

Epidermis



•660nm Rotlicht (Sichtbar)

Dringt in das oberflächliche Hautgewebe ein, verbessert die Durchblutung und reduziert die Wundentzündung.

Dermis




•850nm Nahinfrarotlicht (Nicht Sichtbar)

Dringt tief in das Hautgewebe ein, um tiefe Muskelschmerzen zu lindern und die Gewebereparatur zu fördern.

Hypodermis

Muscle



## Verwendung der Rotlichtlampe



Besserer Schlaf



Verringerung von Gesichtsfalten



Verbesserte Blutzirkulation



Reduziert zelluläre Entzündungen



Verbessert Schuppenflechte



Fördert die Vitalität des Körpers

